

## Description:

GelRed from Biotium is an ultra-sensitive, extremely stable and environmentally safe fluorescent nucleic acid dye designed to replace the highly toxic ethidium bromide (EB) for staining dsDNA, ssDNA or RNA in agarose gels or polyacrylamide gels. It is far more sensitive than EB without requiring a destaining step.

GelRed and EB have virtually the same spectra, so you can directly replace EB with GelRed without changing your existing imaging system. GelRed(TM) can be used to stain dsDNA, ssDNA or RNA in agarose gel via either precast or post gel staining. It can also be used to stain dsDNA, ssDNA or RNA in polyacrylamide gel via post gel staining and it is also compatible with downstream DNA manipulations such as digestion with a restriction enzyme, Southern blotting techniques and clonings.

A series of safety tests have confirmed that GelRed(TM) is noncytotoxic, nonmutagenic and nonhazardous at concentrations well above the working concentrations used in gel staining. As a result, GelRed can be safely disposed in regular trash, providing convenience and reducing cost in waste disposal.

**Concentration:** supplied as 10,000X solution in water.

## Staining Protocols

### 1. Staining DNA by Post Gel Staining

a. Run gels as usual according to your standard protocol.

b. Dilute the GelRed™ 10000X stock solution about 3,300 fold to prepare a 3X staining solution in water with 0.1M NaCl (e.g. add 15µL of GelRed™ stock solution and 5mL NaCl to 45mL water).

While GelRed™ 1X staining solution can also be used for post gel staining, the sensitivity is generally less than with 3X staining solution

(**NOTE:** inspect the 10000X vial carefully. If dye precipitation occurs, heat up or sonicate the vial).

Use of NaCl in the staining solution is optional. NaCl in the solution enhances the staining, but may promote dye precipitation if the staining solution is to be used repeatedly. Any staining solution to be re-used is preferably stored at room temperature in a dark place to reduce possible dye precipitation problems.

c. Carefully place the gel in a suitable container such as a polypropylene container. Gently add sufficient amount of the 3X staining solution to submerge the gel.

d. Agitate gel gently at room temperature for about 30 minutes. Optimal staining time may vary somewhat depending on the thickness of the gel and the percentage of agarose. For polyacrylamide gels containing 3.5 – 10% acrylamide, typical staining time is 30 minutes up to 1 hour with gels of higher acrylamide content requiring longer staining time. The staining solution can be re-used at least 2-3 times.

The unused staining solution can be stored at room temperature in a dark place.

e. View the stained gel with a standard transilluminator (302nm o 312nm) and photograph the gel using an ethidium bromide filter. Similarly, a SYBR(R) or GelStar(R) filter may also be used for photographing with equally good results.

## 2. Staining of DNA by Precasting GelRed™ Gels\*

a. Prepare agarose gel solution using your standard protocol.

b. Dilute GelRed™ 10000X stock solution into the agarose gel solution at 1:10000 (e.g. 5µL stock solution added to 50mL of gel solution). Since GelRed™ is generally thermally stable the 10000X stock solution can be added while the gel solution is still hot. Make sure that the dye is thoroughly mixed with the gel solution by swirling, stirring, or inversion (**NOTE:** inspect the 10000X vial carefully. If dye precipitation occurs, heat up or sonicate the vial).

Alternatively the stock solution may be pre-combined with agarose powder and electrophoresis buffer of your choice followed by microwaving or other heating procedures commonly used for preparing agarose gels. GelRed™ is compatible with all commonly used electrophoresis buffers.

c. Cast the gel and allow it to solidify. Any leftover gel solution may be stored and re-heated later for additional gel casting. Since GelRed™ is hydrolytically stable, GelRed™ precast gels may be prepared in large quantities and stored for later use. To avoid mould formation, we recommend that the precast gels be stored in a refrigerator.

d. Load samples and run the gels using your standard protocol.

**NOTE:** Use only 1/3 or less of your normal DNA-marker amount to avoid overexposure (no bands of marker visibly separated).

e. View the stained gel with a standard transilluminator (302nm o 312nm) and photograph the gel using an ethidium bromide filter. Similarly, a SYBR(R) or GelStar(R) filter may also be used for photographing with equally good results. **NOTE:** If you consistently see band smearing and/or poor band separation, run a gel and post-stain by following the protocol provided (point 1) to confirm if the problem is caused by the dye or other factors unrelated to the dye. If post gel staining is normal and the problem is not caused by the dye, try any of the followings:

- Lower the amount of nucleic acid loaded.
- Lower running voltage
- Lower the amount of agarose in the gel
- Run a longer gel
- Increase the thickness of the gel
- Increase the solidification time to ensure sharp well formation
- Improve your loading technique or select post gel staining

**\* Precasting GelRed™ gel is not suitable for acrylamide gels. Use post gel staining for acrylamide gels.**

## Features of GelRed™

### Safer than EB

Shown by the Ames test and other tests to be nonmutagenic and noncytotoxic

### Easy disposal

Passed environmental safety tests for direct disposal down the drain or in regular trash

### Ultra-sensitive

Much more sensitive than EtBr and SYBR Safe

### Extremely stable

Available in water, stable at room temperature for long-term storage and microwavable

### Simple to use

Very simple procedures for either precast and post gel staining

### Perfectly compatible with a standard UV transilluminator

GelRed replaces EtBr with no optical setting change

### Perfectly compatible with downstream applications

Compatible with downstream DNA manipulations such as digestion with a restriction enzyme, Southern blotting techniques and cloning.

## Storage and Stability

GelRed™ is a very stable dye. We recommend that the 10000X solutions are stored at RT in a dark place. The solution may also be stored at a lower temperature such as 4°C. Dye precipitation may occur during prolonged low temperature storage. When this occurs, heat up the solution in a hot water bath at 45°C to 50°C for two minutes and vortex the solution.

The 1X or 3X post column staining solutions may also be stored at room temperature in a dark place for at least one year. Exposure to light should be avoided during long-term storage. However, the dye can be handled under ambient light without any problem during staining experiments.

## Toxicity

GelRed™ was subjected to a series of tests by three independent testing services to assess the dye's safety for routine handling and disposal.

These tests include:

1. Glove penetration test.
2. Cell membrane permeability and cytotoxicity test.
3. Ames test.
4. Environmental safety tests. Test results confirm that the dye is impenetrable to both latex gloves and cell membranes. The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. GelRed™ appears to be completely cell membrane-impermeable, which may be a key factor responsible for the observed low toxicity. However, since these tests were not performed on human, we still advise that researchers exercise precautions when handling the dye or any other DNA-binding molecules by wearing protective gears.

## Disposal

GelRed™ has successfully passed environmental safety tests in compliance with CCR Title 22 Hazardous Waste Characterization. As a result, GelRed™ is not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash, providing convenience and reducing cost in waste disposal.

**Note:** \*GelRed and its uses are covered by pending US and international patents.

\*\*SYBR is trademark of Molecular Probes, Inc.

## Trouble shooting

### Problem Suggestion

Smeared DNA bands in gel 1.Reduce the amount of DNA loaded by one-half to one-third.

Blown out or smeared bands can be caused by overloading. This is frequently observed with DNA ladders.

2.Perform post-staining instead of pre-casting.

3.Pour a lower percentage agarose gel for better resolution or large fragments.

4.Change the running buffer. TBE buffer has a higher buffering capacity than TAE

Discrepant DNA migration in pre-cast gel

GelRed is designed to be larger than other dyes to prevent it from entering cells, thus rendering the dye safer. The migration of DNA may be affected depending on the dye:DNA ratio.

1.Reduce the amount of DNA loaded by on-half to one-third.

2.Reduce the amount of dye used, i.e. use 0.5X in pre-cast gels.

3.Post-stain the gel in 3X GelRed to avoid any interference the dye may have on migration during electrophoresis.

Weak fluorescence, decreased dye performance over time, or film of dye remains on gel after post-staining. The dye may have precipitated out of solution.

1-Heat GelRed solution to 45-50°C for two minutes and vortex to redissolve.

2.Store dye at room temperature to avoid precipitation.

## Order information

| Prod. No. | Description                    | Quantity   |
|-----------|--------------------------------|------------|
| S9420     | GelRed for Agaros or Page Gels | 0,5 ml     |
| S9425     | GelRed for Agaros or Page Gels | 5 x 0,5 ml |