

Description:

Low Melt Agarose is a specialist agarose for the isolation of nucleic acid fragments between 80 bp and 20 kbp. Low Melt Agarose melts at lower temperatures forming very transparent gels and is suitable for all standard 'in gel' techniques.

Quality:

Molecular Biology Grade

- Certified free of DNases and RNases
- No DNA binding
High lot-to-lot consistency

Order Information

Prod. No.	Description	Quantity
6905-010	Agarose Low Melt	10 g
6905-025	Agarose Low Melt	25 g
6905-100	Agarose Low Melt	100 g