

Description/Preparation:

The Protein Marker US7 is a mixture of 7 purified proteins which are re-dissolved 'ready-to-use' in loading buffer. The proteins resolve into 7 sharp bands in the range of 14.4 kDa to 116.0 kDa when analyzed by SDS-PAGE and stained with optimized dye.

Order Information

Prod. No.	Description	Quantity
3910001	GENE Protein Ladder US7 unstained, 7 bands	2x 1 ml
3910002	GENE Protein Ladder US7 unstained, 7 bands	5x 1 ml