

Description/Preparation:

Protein Marker XXL DeLuxe is a mixture of 10 coloured proteins. These proteins are recombinantly produced in E.coli, highly purified and supplied 'ready to use' in loading buffer. The protein concentrations are optimized to yield well-defined bands directly visible in SDS-polyacrylamide gels. The marker is ready-to-use. There is no need to boil.

Applications:

- Monitoring of protein migration during SDS-polyacrylamide gel electrophoresis
- Monitoring of protein transfer onto membranes during Western blotting
 - Sizing of proteins on SDS-polyacrylamide gels and Western blots

Order Information

Prod. No.	Description	Quantity
3910007	GENE XXL Protein Ladder prestained, Deluxe, 10 bands	2x 250 µ
3910008	GENE XXL Protein Ladder prestained, Deluxe, 10 bands	10x 250 µ