

Synonym	[E.C. 3.2.1.17], Muramidase
State of matter	Solid
origin	from chicken egg white
additional product description	Hydrochloride form salt-free, albumin-free, freeze-dried
Mass	approx. 14400 g/mol
CAS-No.:	9001-63-2
EC-No.:	232-620-4
Storage:	-20°C
Shipment:	RT
LGK:	10 - 13
WGK:	1
Specification:	
DNases/RNases/Proteases	not detectable
Activity (pH 6.2)	min. 20000 U/mg

Order Information

Prod. No.	Description	Quantity
CE9188	Lysozyme	1 g
CE9189	Lysozyme	10 g