

Description:

Maximo Dry Master-Mix is optimized and ready-to-use mixture of all components for a successful PCR. Only diluent-solution your primers and DNA Template has to be added.

Maximo Dry Master-Mix contains a thermostable DNA polymerase that possesses a 5' → 3' polymerase activity and a double-stranded specific 5' → 3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94kDa.

Features:

- All components: **"Hot-Start"** Taq polymerase, dNTP's, reaction buffers, a blue tracking dye, enhancers and stabilizers are lyophilized
- Transportation and short time storage at room-temperature
- Simple, and fast setting up procedure for high yield and repeatable results
- Direct-loading on the Agarose-gel saves time

Order Information

Prod. No.	Description	Quantity
S9295	Maximo Dry Master Mix	12x 8 flat-cap