

Description:

Maximo Taq DNA Polymerase is a thermostable DNA polymerase that possesses a 5' → 3' polymerase activity and a double-strand specific 5' → 3' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of 94kDa.

Features:

Maximo Taq DNA Polymerase provides robust PCR performance in a wide range of PCR applications and different templates. Best value in terms of cost per unit.

Order Information

Prod. No.	Description	Quantity
S9101	Maximo Taq DNA Polymerase	500 Units
S9102	Maximo Taq DNA Polymerase	5x 500 Units
S9103	Maximo Taq DNA Polymerase	20x 500 Units