

PCR Detection-Test-systems for the identification of infection diseases provide a convenient analysis of infection agents presence in biological samples.

All reagents (Taq DNA polymerase, dNTPs, specific primers, salts and stabilizers) are lyophilized in PCR tubes (0,2 ml), only sample to be analyzed and dilution buffer should be added before running the PCR.

The PCR-Test-System contains:

- 80 DNA test tubes
- 10 positive controls
- 10 negative controls
- in total: 100 tubes with lyophilized reagents (1 tube – 1 reaction)
- Specific dilution buffer

Order Information

Prod. No.	Description	Quantity
M92115	Mycobactericum	80 + 10 + 15
M92016	Mycobactericum avium, subsp.	80 + 10 + 16
M92026	Mycobactericum intracellulare	80 + 10 + 17
M92015	Microbactericum tuberculosis	80 + 10 + 18