
Synonym	β -Nicotinamide Adenine Dinucleotide
---------	--

State of matter	Solid
Formula	$C_{21}H_{26}N_7Na_4O_{17}P_3$
M	833.35 g/mol
CAS-No.:	2646-71-1
EC-No.:	220-163-3
Storage:	-20°C, protected from light
LGK:	10 - 13
WGK:	1
Specification	
Assay (enzym.)	min. 95 %
Water (K.F.)	max. 8.0 %
UV Absorption:	
A340/A260 (pH 10)	0.43 \pm 0.01

Order Information

Prod. No.	Description	Quantity
CE9202	NADPH - Tetrasodium salt	25 mg
CE9203	NADPH - Tetrasodium salt	100 mg
CE9204	NADPH - Tetrasodium salt	500 mg