

## Ensure High Yield Protein Expression

Transient transfection of suspension-growing cells represents an attractive alternative to costly, time-consuming and labour intensive stable transfection.

The ProteinTrans transfection Kit has been designed for large scale up transient transfection and allows high protein production yield with reproducible results.

The ProteinTrans transfection Kit is dedicated to achieve High Yield Protein Expression in mammalian cells. This Kit has been designed for maximum efficiency in HEK293 and CHO cells growing in suspension.

The ProteinTrans transfection Kit, is ideal for bioreactor, spinner or flasks.

- **High yield recombinant protein production**
- **Suitable for both HEK293 and CHO cells growing in suspension**
- **Compatible with any synthetic or regular media used for protein production**
- **Completely animal origin free**

## Order Information

Prod. No.	Description	Quantity
SB901500	Protein Trans Reagent + CHO-Blast	1.5 ml + 5 ml
SB903000	Protein Trans Reagent + CHO-Blast	3 ml + 10 ml
SB915000	Protein Trans Reagent + CHO-Blast	15 ml + 50 ml
SB930000	Protein Trans Reagent + CHO-Blast	30 ml + 100 ml