

<b>Synonym</b>	<b>Carbamide, Carbonyldiamide</b>
State of matter	Solid
Solubility (20°C)	1000 g/L (H <sub>2</sub> O)
Melting point	133°C
Formula	NH <sub>2</sub> CONH <sub>2</sub>
Mass	60.06 g/mol
CAS-No.:	57-13-6
EC-No.:	200-315-5
Storage:	RT
LGK:	10 - 13
WGK:	1
<b>Specification:</b>	
DNases/RNases/Proteases	not detectable
Assay (from N)	min. 99.5 %
pH (8 M; H <sub>2</sub> O; 20°C)	7.5 - 10.5
Heavy metals (as Pb)	max. 0.001 %
Chloride	max. 0.0005 %
Sulfate	max. 0.001 %
Fe	max. 0.0005 %

## Order Information

<b>Prod. No.</b>	<b>Description</b>	<b>Quantity</b>
CE9167	Urea Crystalline	1000 g
CE9168	Urea Crystalline	5000 g