

Specific transfection reagent for Vero Cells

VeroTrans is a powerful reagent dedicated to the transfection of Vero cells.

VeroTrans is specifically designed to obtain highly efficient and reproducible transfection of Vero Cells. It can be used for many applications such as stable and transient transfection, protein and viral production.

High transfection efficiency in Vero cells

VeroTrans reagent has an exceptional capacity in destabilizing cell membranes, allowing the delivery of important DNA amounts into cells. Transfection efficiency exceeds most of others transfection reagents.

Order Information

Prod. No.	Description	Quantity
S960250	VeroTrans	250 µg
S960500	VeroTrans	500 µg
S961000	VeroTrans	1 ml
S965000	VeroTrans	5x 1 ml