

Description/Preparation:

The pBR322 DNA/Hinf I digest yields the 10 fragments. pBR322 DNA was completely digested by Hinf I, phenol/chloroform extracted, ethanol precipitated, dissolved in 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.

Features:

Supplied in 6X Loading buffer

Order Information

Prod. No.	Description	Quantity
3900040	pBR322 DNA / Hinf I digest	200 µ
3900041	pBR322 DNA / Hinf I digest	5x 200 µ